

UNICARE (GS003*)

Disposable-Examination gloves made of nitril
 powder free • blue • latex free • Chemical protective
 glove EN374 • non-steril • suitable for short-term food contact

Specifications:

Material: Nitril
 Powder: powder free
 Cuff: beaded cuff
 Design: ambidextrous
 Surface: Microroughened
 Farbe: Blue
 Quality control:
 PSA – EN 374 CAT III - TypC,
 EN 388, EN 420; MPG – EN
 455 Teil 1–4
 EG-Standard 596/2009
 EG-Standard 202/2014; EN
 ISO 9001:2008 EN ISO
 13485:2012

Item Number and Size:

GS0031 | Extra Small
 GS0032 | Small
 GS0033 | Medium
 GS0034 | Large
 GS0035 | Extra Large

Packaging:

100 pieces per Box
 10 Boxen per Carton

Classification:



For more
 Information
 click on
 the
 image



Test parameter

Test criteria	Level	AQL
Dimensions	S-2	4,0
Physical Properties	S-2	4,0
1000-ml-Watertightness test	G-1	1,5
Visual inspection for major errors	G-1	2,5

Dimensions

Item.-No.	GS0031	GS0032	GS0033	GS0034	GS0035
Size	XS = 5–6	S = 6–7	M = 7–8	L = 8–9	XL = 9–10
Width	75 ± 5 mm	80 ± 10 mm	95 ± 10 mm	110 ± 10 mm	≥ 110 mm
Length	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm	≥ 240 mm

Material thickness

Measuring	One layer (mm)
Cuff (at a distance of 40mm to the beaded edge)	0.06 ± 0.02
Palm (in the middle of the palm)	0.10 ± 0.03
Fingertips (at a distance of 15mm to the tip of the middle finger)	0.11 ± 0.03

Physical Properties

Parameter	Before aging	After aging
tensile strength (N)	min. 6.0	min. 6.0
Elongation at break (%)	min. 500	min. 400
modulus of elasticit (MPa)	min. 14.0	min. 14.0

Chemical protection

The product complies with the PSA regulation (EU) 2016/425

Type-Examination (Modul B) & Modul C2

SATRA Technology Europe, Bracetown Business Park, Clonee, D15 YN2P Ireland
 (Notified Body 2777)

Test chemical	EN ISO 374-1:2016 Permeationslevel	EN374-4:2013 Degradation %
Natriumhydroxid 40% (K)	6	-78.9

Hint

The above information applies to the standard version.
 Subject to changes

