


SMART PRO - ADVANCED CONTROLLER 7" COLOUR TOUCH PANEL GUARANTEES INTUITIVE AND COMFORTABLE OPERATION



ADVANTAGES OF THE SMART PRO CONTROLLER

- large (7"), clear, full-colour touch screen
- LAN, USB ports, Wi-Fi for communication and data transfer
- multi-segment time and temperature programs
- overview of data in tabular and graphic form
- visual and audible alarm
- administration functions for easy management
- password protected log-in
- internal memory for programs and data storage
- operating with gloves on
- event registry with user notifications
- LabDesk software and instruction manual for direct download
- Alarm Bar – instant visual information about chamber status
- Quick Note – user can save text notes (50 characters) in Smart PRO controller memory
- Quick Change of parameters: temperature, humidity, time, air flap and fan (according to model)